



CHANGING YOUR CPAP/BiPAP SUPPLIES

TIP: For this schedule, "Month #1" should be the CURRENT month. "Month #2" is NEXT month, etc...

Replacement		Schedule											
Item	Frequency	MONTH 1	2	3	4	5	6	7	8	9	10	11	12
Mask													
___ Full Face Mask 	1 every 3 months	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>
___ Nasal Mask 		--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>
___ Nasal Pillow Mask 		--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>
Mask Cushions & Pillows													
___ Full Face Mask Cushion 	1 per month	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
___ Nasal Mask Cushion 	2 per month	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
___ Nasal Pillows 		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Headgear	1 every 6 months	--	--	--	--	--	<input type="checkbox"/>	--	--	--	--	--	<input type="checkbox"/>
Chinstrap 	1 every 6 months	--	--	--	--	--	<input type="checkbox"/>	--	--	--	--	--	<input type="checkbox"/>
Tubing 	1 every 3 months	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>	--	--	<input type="checkbox"/>
Disposable Filters 	2 per month	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Non-disposable Filters 	1 every 6 months	--	--	--	--	--	<input type="checkbox"/>	--	--	--	--	--	<input type="checkbox"/>
Humidifier Water Chamber 	1 every 6 months	--	--	--	--	--	<input type="checkbox"/>	--	--	--	--	--	<input type="checkbox"/>

*Quantity and frequency may vary with your insurance coverage.

Sleep Central is a division of Rotech Healthcare Inc. Sleep Central specializes in providing the highest quality supplies for your CPAP/BiPAP and the ongoing customer support needed for your success and improved health. Respiratory Therapists are available to answer technical questions or clinical issues and Patient Care Coordinators are available to address your questions and monitor your replacement schedule inline with your insurance coverage.

Sleep Central: (844) 727-6275 or email help@sleepcentral.com