

Stop Worrying about Running Out of Oxygen.



Oxy
GO
Keep Going.™



www.rotech.com

1-866-335-2955

What is OxyGo?

The OxyGo® family of portable oxygen concentrators (POC) is designed to keep oxygen patients—like you—GOING!

An OxyGo POC is a lightweight medical device that takes the air around you and converts it into the oxygen prescribed by your doctor.

The OxyGo family has two members: **OxyGo** and **OxyGo FIT™**

Both are industry leaders in reliability, usability and quality. And both are clinically validated to be used 24/7... so they can supply the medical oxygen you need at home or on the go.

OxyGo is different from other POCs because it shows you exactly how much time you have left on the battery, and you can change the battery easily. You can also charge OxyGo and use at the same time—either in the car or at home!

Will OxyGo work for me?

Depending on your oxygen prescription, it's likely that it will! The most common prescription for oxygen is 2 liters per minute, and either OxyGo or OxyGo FIT can deliver that amount of oxygen—and *more!*—in pulse dose.

The OxyGo and the OxyGo FIT are available by prescription only. If you have been prescribed oxygen, there is a good chance that either member of the OxyGo family will work for you. However, you will need your oxygen provider to perform an oxygen assessment to see if you can saturate with OxyGo.

How do I get OxyGo?

Renting or owning your OxyGo, OxyGo FIT and accessories may be easier than you think! Many insurance plans will cover the rental of the OxyGo. However, you may choose to keep your OxyGo forever by using cash, credit or debit cards... or one of our available, convenient financing options.

Call 1-866-335-2955 to get started!



OxyGo shows you how much time is left before you need to recharge or switch batteries

OxyGo® is ready to go when you are!

OxyGo keeps you active. Whether you are out on the town or home with family — you can always have your oxygen with you. Learn more at: www.oxygo.life



OxyGo FIT™

5.9" x 2.7" x 7.2" 2.8 lbs.

We named our newest POC the OxyGo FIT because it is going to "fit" into your lifestyle! It weighs just 2.8 pounds and has 3 settings, with 210 ml of oxygen per setting. The OxyGo FIT is well constructed and will operate for almost 3 hours using a standard, user-replaceable battery. It also has an optional larger battery, which will power the unit for up to 5 hours* and adds less than one-half pound in weight.

OxyGo FIT (with 4-cell battery)

#1400-2000 MAP = \$2,495

- 2.8 lbs. | 5.9" x 2.7" x 7.2"
- 3 Pulse Flow Settings (630 ml/min)
- Up to 5 hour battery life*
- Only 40 dB. at setting 2
- Included 3-year warranty on the unit, 1-year warranty on sieve beds & battery

* double cell battery, setting 1

OxyGo FIT comes with:



- Standard, 4-cell battery
- DC car charger, and AC wall charger
- Nasal cannula
- Carry bag with padded, adjustable strap

Optional Accessories



Additional Batteries

8-Cell	4-Cell
#1400-2010-8	#1400-2010-4
MAP = \$495.00	MAP = \$325.00



Desktop Battery Charger

#1400-2030
MAP = \$275.00



Air Cabinet Filter

(1 per package)
#1400-2310
MAP = \$11.95



Backpack

#1170-2420
MAP = \$100.00



OxyGo 5 Setting

8.75" x 3.0" x 7.25" 4.8 lbs.

The original OxyGo has 5 pulse flow settings (highest portable O₂ output available today), 1,050 ml output, and up to 10.25 hours* of battery life—all in a compact, portable unit that weighs less than 5 pounds. OxyGo also has user-replaceable batteries and sieve beds, so it's ready when you need it. And with its NEW 5 year warranty, OxyGo is truly designed to help you Keep Going.

OxyGo 5 Setting (with 8-cell battery)

#1400-1000 MAP = \$2,495

- 4.8 lbs. | 8.75" x 3.0" x 7.25"
- 5 Pulse Flow Settings (1050 ml/min)
- Up to 10.25 hour battery life*
- Only 39 dB. at setting 2
- Included 5-year warranty on the unit, 1-year warranty on sieve beds & battery

* double cell battery, setting 1

OxyGo comes with:



- Standard, 8-cell battery
- DC car charger, and AC wall charger
- Nasal cannula
- Carry bag with padded, adjustable strap

Optional Accessories



Additional Batteries

16-Cell	8-Cell
#1400-1010-16	#1400-1010-8
MAP = \$495.00	MAP = \$325.00



Desktop Battery Charger

#1400-1030
MAP = \$275.00



Air Cabinet Filters

(2 per package)
#1400-1311
MAP = \$29.00



Backpack

#1170-1420
MAP = \$100.00



Get the medical equipment you need **TODAY.**

We are Rotech Healthcare—a family of hometown healthcare companies and industry leaders in providing high quality medical products, unsurpassed customer service, and assistance in patient home health treatment. With hundreds of locations from coast to coast, you will find a local Rotech company near you—staffed with good neighbors... good friends...and a team who genuinely cares about patients.

Each and every day, thousands of dedicated Rotech employees bring our special brand of care and quality service to neighborhoods and communities across the U.S. We offer a comprehensive selection of home medical equipment and services to help patients better manage their treatment and care at home.

Our constant goal is to enable you to manage your care and live a more comfortable and productive life—every day.

CALL NOW

**And let our dedicated team help you
get the equipment and care you need!**

1-866-335-2955